



Analytical Balance & Scientific Instruments



LMIF SERIES



ICE FLAKER

SALIENT FEATURES

- Fully Automatic microprocessor control
- Compact Design and low maintenance
- Stainless Steel construction
- Continuous Ice Flakes Output
- CFC free compressor
- Noiseless operation
- Overload protection
- Low water level detection
- Wheels for easy mobility

APPLICATION :

Research Institute, Hospitals, Laboratories
Biotechnology & Universities.

SPECIFICATIONS

MODEL	LMIF-30	LMIF-50	LMIF-85	LMIF-100
Make	LABMAN			
Capacity of Ice flake (kg per 24 hrs.)	30 kg	50 kg	85 kg	100 kg
Storage bin	ABS Plastic			
Storage Capacity (kg)	10 kg	15 kg	25 kg	
Cooling Method	Air			
Working temperature	Room temp should be below 30 °C			
Construction	Stainless Steel			
Noise level	52 dB			
Water Supply	3 / 4 inch BSP			
Refrigerant	R134A			
Input Power	180W	280W	420W	
Dimension (mm)	300X493X547	380X543X722	548X611X883	
Power supply	220V / 50Hz			
Gross Weight	36kg	45kg	62kg	64kg