



Analytical Balance & Scientific Instruments



LMWB SERIES

SALIENT FEATURES



DIGITAL WATER BATH

Novel double walled structure effectively prevents vapor leakage

Intelligent temperature controller with PID program

Digital display and touch operating buttons

Function of timing and over temperature protection

Stainless Steel inner chamber .

Applications :

The equipment applies to Medical Unit, college, Scientific Research Institution and The Department of Chemical Printing, Dyeing, Pharmacy, Petroleum and Electronic Component Production for the use of evaporation, concentration and constant temperature heating.

SPECIFICATIONS

MODEL	LMWB-02	LMWB-04	LMWB-06	LMWB-08
Make	LABMAN			
Temperature Uniformity	± 0.5			
Temperature Range	Ambient +5°C-100°C			
Inner Chamber dim. (SS) WxDxH (mm)	168 x 168 x 120	325 x 325 x 120	480 x 325 x 120	635 x 325 x 120
Outer dimension (mm)	360 x 240 x 190	510 x 380 x 190	670x 380 x 190	830 x 380 x 190
Capacity (approx.ltr.)	3.5	13	19	25
Power Wattage	400 W	1000 W	1500 W	2000 W
Power supply	220 V / 50 Hz			
Gross Weight	4.5 kg	8 kg	9.5 kg	10 kg