



## Analytical Balance & Scientific Instruments



### LMGG Series



### HYDROGEN GAS GENERATOR

### SALIENT FEATURES

- High Purity Hydrogen Gas Generator
- Hydrogen is produced by the electrolysis of water
- The electrolysis separation pool is made of S.S
- Low temp. of electrolysis pool.
- Equipped to prevent the returning of the liquid.
- Flow is visualized by digital display
- Over voltage & over current protection

### SPECIFICATIONS

MODEL	LMGG-300H <sub>2</sub>	LMGG-500H <sub>2</sub>	LMGG-5000A
<b>TYPE OF GENERATOR</b>	<b>HYDROGEN GAS GENERATOR</b>		<b>AIR GENERATOR</b>
Make	LABMAN		
Purity	99.999%		
Output Flow(H <sub>2</sub> ) ( Max.)	300ml/min	500ml/min	N/A
Output Flow (Air)	N/A	N/A	5L /min
Output Pressure	0.4Mpa		0.6Mpa
Ambient temperature	0 – 40°C		
Power Wattage	150 W		
Gen. Dim. LxWxH (mm)	350 x 190 x 300		
Power supply	220 V / 50 Hz		
Gross Weight	14 Kg		15 kg