



Analytical Balance & Scientific Instruments



SALIENT FEATURES

- Sturdy & Robust Stand
- Anti Rust material.
- Dial Diaphragm
- Coarse Focusing Adjustable Mechanism
- Built in LED Illumination

MONOCULAR MICROSCOPE

SPECIFICATIONS

MODEL	LMSC-38M
Make	LABMAN
Objective	Achromatic 4 X, 10X, 40X(S)
Viewing Bodies	Fixed Monocular Head 38° Inclined
Eyepiece	Wide field Plan scope eyepiece WF 10X
Nose Piece	Six stop Dial Diaphragm
Mechanical Stage	Plain stage with paired clips 90 x 90mm
Condenser	Single Lens N.A. 0.65
Focusing	Coarse Focusing Adjustable 8mm
Illumination	LED
Power supply	220V/50Hz
Gross Weight	2.5 kg