



Analytical Balance & Scientific Instruments



LMIO SERIES



INCUBATOR

SALIENT FEATURES

- High-quality cold rolling steel electrostatic spraying exterior.
- Stainless steel inner chamber, foursquare semicircle transition convenient for cleaning.
- Adjustable shelves
- Micro-computer intelligent control system with PID program.
- Large screen digital display and soft touch operation keys.
- Over temperature protection and auto-power break-off.
- Digital Timer
- Double door seal.

SPECIFICATIONS

Model	LMIO-43	LMIO-124	LMIO-226
Make	LABMAN		
Temperature Range	Ambient + 5 to 65 °C		
Temperature Resolution	0.1°C		
Temperature Uniformity	± 0.5 °C		
Working temperature	5°C to 40°C		
Heating Power	250 W	500 W	600 W
Inner Chamber Dimension (SS) WxDxH (mm)	350x350x350	500x450x550	600x580x600
Outer Dimension (mm) WxDxH	485x480x570	626x590x755	726x720x810
Digital timer	1-9999 min		
Load per rack	2 /15 kg per rack		
Power Supply	220V/50HZ		
Gross Weight	34 kg	46 kg	54 kg