



Analytical Balance & Scientific Instruments



SALIENT FEATURES

- Stainless steel inner chamber;
- Four square semicircle transition, convenient to clean;
- Micro-computer control system; LCD display with functions of time setting, parameter memorizing and power reset recovering.
- CFC-free refrigeration and gentle breeze circulation.
- Double door seal, tempered glass, silicon rubber seal, Inner & Outer door with magnetic rubber seal.
- Inner chamber with light for convenient observation.

BOD INCUBATOR

SPECIFICATIONS

MODEL	LMBI-150
Make	LABMAN
Temperature Range	0 - 65°C
Temperature Resolution	0.1°C
Temperature Uniformity	± 0.5°C
Working temperature	5°C to 30°C
Timer	1-9999 min
Refrigerant	R134A
Power Wattage	500 W
Inner dimension (mm) WxDxH	500 x 400 x 750
Outer dimension (mm) WxDxH	640 x 620 x 1270
Load per rack	2 / 15kg
Power Supply	220 V / 50 Hz
Gross Weight	97 kg