



Analytical Balance & Scientific Instruments



LMVO Series



VACUUM OVEN

SALIENT FEATURES

Stainless steel Inner chamber & Outer body electro static spraying exterior

Micro-computer temperature controller with function of timing and over-temperature protection.

High brightness digital tube display and soft touch keys for parameter settings.

Tightness of door is adjustable

Double layered tempered glass anti-explosion door; silicon rubber seal provides a high vacuum degree and keep longer time.

SPECIFICATIONS

Model	LMVO-24	LMVO-52
Make	LABMAN	
Temperature Range °C	Ambient + 10- 250°C	
Temperature Resolution °C	0.1°C	
Temperature Uniformity	± 1%	
Inner Dimension (mm) WxDxH	300x300x270	415x370x340
Outer Dimension (mm) WxDxH	480x480x605	560x540x680
Power Wattage	800 W	1400 W
Vacuum degree	133pa	
Power supply	220V / 50Hz	
Gross Weight	52 kg	92 kg